

Canine Rehab Systems

Specifications Sheet

Mobility Frame

Model MF-01

GENERAL:

This frame is designed to be an integral part of a rehabilitation practice that works with mobility impaired dogs. It allows the therapist to safely and comfortably support the dogs with assisted ambulation activities.

The open frame allows easy access to all parts of the dog's body and extremities, allowing the therapist to easily provide neuromuscular stimulation or other techniques for the patient. The caster wheels lock, making the transfer in and out of the frame more secure.

Five different slings (XS-S-M-L-XL) allow for secure support of dogs of all. Each sling has adjustable straps, to allow for dogs of different height to correctly fit into the frame. As the dogs get stronger, the straps can be lengthened, allowing the dogs to bear more weight through their limbs, but still preventing them from falling or collapsing. These slings are also interchangeable with our TreadRails System, used to support the dog while in the underwater treadmill.

DESIGN FEATURES:

The frame is made of a strong extruded aluminum, with adjustable eye bolts which allow accurate positioning of the dog in the frame. The slings are made of a durable yet soft nylon fabric which can be used both on land, and in the water. The frame and harness will support dogs up to 250 lbs., so it can be used with dogs of all sizes. The modular frame can be easily disassembled for storage or transport, and reassembled in minutes. Each Canine Standing Frame is custom manufactured - dimensions may be specified to fit individual clinic needs.



Canine Rehab Systems LLC

CanineRehabSystems.Com



3787 Jacks Creek Rd NW

Monroe, GA 30655

Info@CanineRehabSystems.Com